


Cyanoacrylate Treated Fingerprints

Cyanoacrylate treated fingerprints dyed with RAM on metal foil:



No Light, No filter



BMT0 Head with amber filter

Cyanoacrylate treated fingerprints dyed with RAM on plastic:




No Light, No Filter



BMT0 Head with amber filter

Cyanoacrylate treated fingerprints dyed with Rhodamine 6G on plastic:



No Light, No Filter




BMT0 Head with amber filter

Cyanoacrylate treated fingerprints dyed with Basic Yellow on plastic:




No Light, No Filter




BMT^o Head with amber filter

Cyanoacrylate treated fingerprints dyed with MBD on plastic:



No Light, No Filter



BMT^o Head with amber filter